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Manufactured in P.R.C., for FEMIEKO APS, Ryttervangen 11, 7323 Give, Denmark.

Distributed with love by FEMIEKO GROUP S.R.L., DV24, CI-24, Bucharest, Romania.

DAGH Distribution: Abraxa Business Network GmbH, Am Schafberg 2a. 96489 Niederfullbach, Deutschland.

## MODERN SOLUTION FOR INTIMATE BALANCE

### ABOUT UNIVERSAL SAFE MENSTRUAL DISC

The **Femi.Eko® Menstrual Disc** is a flexible, disc-shaped device that is placed at the base of the cervix to collect menstrual blood. It offers a discreet and modern menstrual experience and is a comfortable and safe alternative to conventional menstrual products. It is made from 100% biocompatible, medical grade Elastosil silicone from Wacker Labs Germany and contains Saniconcentate™, an innovative medical grade compound patented by Parx Materials. What makes Saniconcentate™ truly remarkable is its ability to provide 99.99% natural protection against bacteria, microbes and viruses without migrating or dissolving during use, and is commonly used in implantable medical devices. This means you do not have to worry about harmful chemicals, bacteria, irritation and therefore the need for frequent sterilisation. This is because the antimicrobial properties are activated by simply rinsing under water, so the disc is safe to use immediately after purchase. Specially designed to fit your anatomy, this disc moulds perfectly to the vaginal canal and protects you continuously for up to 12 hours and can be reused for at least 5 years.

The **Safe Universal Menstrual Disc with Saniconcentate™** is a Femi.Eko® innovation approved in Europe which, in addition to its antibacterial properties, has an integrated removal ring that facilitates removal and improves the experience of using this menstrual disc.

### HOW DOES IT WORK?

The Universal Safe Menstrual Disc is designed to fit comfortably in the vagina and is taking an anatomical shape. It collects the menstrual blood inside, under the cervix, i.e. it's placed as deep and towards the coccyx as possible. The disc is located at the end of the vaginal canal, just below the cervix. It's therefore very important to push it in as far as possible and keep it in place by positioning the outer edge of the disc, the one with the ring, behind the bone at the front of the pelvis, known as the 'pubic bone'.

When inserted correctly, the disc moulds so comfortably into the vagina that you won't feel it at all and you'll even forget you're your period. By correct positioning, we mean that the disc is placed at an oblique angle inwards so that it's in the right place for optimal function and comfort.

Unlike menstrual cups, it doesn't rely on the creation of a vacuum to stay in place, making it much easier and more user-friendly, but especially suitable for IUD users.

### HOW TO USE?


The **Femi.Eko® menstrual disc** does not need to be sterilised after purchase, but we still recommend that you sterilise it before first use to be on the safe side :). You can sterilise the disc by boiling it in boiling water for 3 to 5 minutes (make sure the disc does not touch the bottom of the boiling pot) or by using a microwave, UV or other steriliser from [www.emieko.com](http://www.emieko.com)

After this initial sterilisation, further sterilisation between uses is not necessary. Simply rinse it under a stream of water to activate the Saniconcentate™ antibacterial compound.

To use the menstrual disc successfully, we recommend that you follow the steps below:


Hygiene is always important when it comes to the genital area and not just there! So wash your hands thoroughly before using the disc. Also wash the Safe menstrual disc with water and/or oil-free, pH-neutral soap. Rinse it well!


Fold the Femi.Eko® menstrual disc into a figure 8 and position the extraction ring towards the outside of the vagina. Can you see it? It's no bigger than a tampon, so relax, it'll be a piece of cake :) If you wish, you can

use water or a water-based lubricant (not silicone) to make insertion easier.


Hold the disc firmly between your thumb and forefinger and start pushing it slowly through the opening of the vulva. As soon as you have passed the vaginal entrance, push the disc further (still holding it firmly in your hand) and tilt it towards the coccyx. Push as deeply as possible until you feel that the disc has

nowhere to go.

Once you have pushed the disc in deeply, it's important to fix it above the pubic bone. Use your index finger to lift the edge of the ring disc towards your navel. Once you have carried out these steps, you can forget about menstruation for 12 hours (many women even report fewer cramps with the disc).

To remove the disc, you must be relaxed. So don't forget to breathe :) The disc is there and it just takes a little skill. Insert your finger into your vagina and grasp the special ring of the disc or the thickened rim. You can do this with either finger, whichever is more comfortable for you. Grasp the ring and carefully

pull the disc out. If you cannot reach the disc, read the following suggestions.

Empty the disc (we recommend that you do this immediately after flushing it down the toilet by simply grasping the ring without running the risk of dropping it in the toilet), rinse it well with cold water and then wash it

with neutral soap and warm water, rinse it very well. Then you can use it again. When the cycle is finished, wash the disc and store it in the special **Femi.Eko® bag** you received in the package.

### SUGGESTIONS FOR USE:

- ✿ Empty it at least twice a day (every 12 hours); we recommend before going to bed and after waking up.
- ✿ Always wash your hands and menstrual disc thoroughly before and after use.
- ✿ Do not share the menstrual disc with other people.
- ✿ Replace the menstrual disc if it is torn or damaged.
- ✿ The menstrual disc is not recommended for use during

bleeding after childbirth or after surgery.

- ✿ You can use the Safe Disc with an intrauterine device or a contraceptive ring, but if in doubt, ask your doctor first.
- ✿ It can be used safely with most contraceptives, but ask your doctor if you have any questions.
- ✿ If you drop your menstrual disc in the toilet, you do not need to replace it. Before using it again, disinfect it by boiling it in water for 4-5 minutes (but no longer than 7 minutes) or wiping it thoroughly with a soft cloth and 70% isopropyl alcohol and allowing the alcohol to evaporate.
- ✿ The menstrual disc does not protect against pregnancy or sexually transmitted infections (STIs).
- ✿ Wash the menstrual disc after sexual intercourse.
- ✿ If the disc moves upwards during arousal or intercourse, allow it to return to a lower position before removing it.
- ✿ If you feel uncomfortable during intercourse, we recommend that you change the position of the menstrual disc or try a different position that is more comfortable for you and your partner.
- ✿ If you cannot get the disc out, squat down and push with your pelvic muscles as if you were having a bowel movement. This movement will change the position of the disc and you'll be able to grasp it by the ring. If you cannot grasp the ring, turn the disc with your finger in the vagina until it has reached a comfortable position. If you cannot grip the ring, grasp the edge of the disc with your fingers and pull gently.
- ✿ Don't hurt yourself, be gentle with your movements and don't insert any objects that could injure your mucous membrane.
- ✿ If you notice any leaks when using the disc, check whether you have inserted it deep enough or whether you have an anatomical feature, such as a retroverted uterus, that requires special adjustment. Make sure that the edge of the disc is positioned correctly over the pubic bone, with the ring facing upwards, to prevent it from slipping (the disc may not sit properly and slip).
- ✿ If you leak into the toilet while urinating or defecating, do not worry. This phenomenon is called auto-dumping and is a feature that many women appreciate. Refit the disc behind the pubic bone after this experience, but remember to replace it at the recommended 12-hour intervals. This phenomenon does not occur during walking or everyday activities, but is a natural process caused by the specific contractions that occur during urination.
- ✿ We recommend that you remove the disc over the toilet or in the shower to avoid splashing blood. If you show signs of menstrual anxiety - menophobia - we recommend that you insert the menstrual disc using the **Femi.Eko® universal applicator** to make the experience easier and reduce contact with the menstrual fluid.

### CARE:

- ✿ Wash thoroughly after each use with warm water and a gentle, pH-balanced intimate soap that is free of silicone, oils, fragrances or chemical ingredients. Avoid harsh cleansers and anything that could irritate the skin. Rinse thoroughly with warm water!
- ✿ Some discoloration may occur over time. This is normal and can vary from person to person. To avoid this, always rinse the disc of blood with cold water first.
- ✿ As a vegan and natural alternative, you can mix a small amount of baking soda or apple cider vinegar with water and

clean the disc with it. Rinse well afterwards and do not leave the disc to soak in these substances.

- ✿ Use a soft brush for dried blood deposits.
- ✿ Occasionally wash and iron the storage bag to ensure complete hygiene.
- ✿ Dry the menstrual disc thoroughly and store it in the cotton bag (or another permeable container). Never use a plastic bag or an airtight container for storage.
- ✿ If the menstrual disc cracks or shows signs of silicone damage (which can occur after using non-recommended cleaning products), it is time to replace the disc.
- ✿ Do not wash the menstrual disc with substances that can irritate the vaginal area. Avoid using cleansers, soaps with strong fragrances or dyes that can disrupt the natural pH or damage the silicone, such as silicone-based products, oils, etc.
- ✿ Do not soak the menstrual disc in any liquid other than water. Silicone is a porous material and can absorb oils or detergents if used over a long period of time. These substances can irritate the vaginal area or damage the silicone.
- ✿ If you want to clean your disc in a public restroom, take a bottle of water with you to clean it and wash your hands thoroughly beforehand.

### SAFETY AND PRECAUTIONS:

Keep the menstrual disc away from children and pets.

Never use the menstrual disc for more than 12 hours at a time.

If you experience symptoms such as general burning, irritation or swelling in the genital area or discomfort when urinating, stop using the disc and consult a doctor.

The menstrual disc does not protect against pregnancy or sexually transmitted diseases.

Internal menstrual products have been associated with toxic shock syndrome (TSS). TSS is a rare but serious condition that can progress rapidly and lead to death. If you experience fever, severe abdominal pain, rash or vomiting while using the disc, stop use immediately and seek medical attention. If you are reading this information, assume that you are aware of the signs of TSS and contact your doctor as soon as possible if anything unusual occurs.

### QUESTIONS?

If you have any further questions or concerns, please do not hesitate to contact us. We are here to assist you and give you the information you need to have the most pleasant and comfortable experience possible. Contact us at: [hello@femieko.com](mailto:hello@femieko.com). We will be happy to answer any questions you may have about the Safe Menstrual Disc.

Thank you for choosing Femi.Eko® for a healthier and greener menstrual experience!

Manufactured in P.R.C., for FEMIEKO APS, Ryttervangen 11, 7323 Give, Denmark.

Distributed with love by FEMIEKO GROUP S.R.L., DV24, CI-24, Bucharest, Romania.

DAGH Distribution: Abraxa Business Network GmbH, Arn Schafberg 2a, 96489 Niederfullbach, Deutschland.

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