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DANISH INTIMATE PRODUCTS
MADE WITH CARE FOR YOU
AND THE PLANET

ABOUT THE MENSTRUAL DISC

A menstrual disc is a flexible, disc-shaped device placed at the base of the cervix to collect menstrual blood and offers a discreet, high-capacity and safe alternative to traditional menstrual products.

Femi.Eko®'s menstrual disc is made from 100% medical grade silicone Elastosil® from Wacker Labs Germany, which is biocompatible and widely used for internal medical devices. This means you don't have to worry about harmful chemicals, bacteria and irritation, promoting intimate balance.

Specially designed to fit every woman's anatomy, this disc is designed to mold perfectly to your intimate area and protect you for up to 12 hours and over 5 years.

The menstrual disc is a unique size and has an integrated removal ring that facilitates and improves the experience of using a menstrual disc, our special design being appreciated by women around the world. Thank you!

HOW DOES IT WORK?

The Femi.Eko® menstrual disc is designed to fit comfortably in the vagina and has a round anatomical shape. It collects the menstrual blood inside, under the cervix, i.e. it is positioned as deep as possible and towards the coccyx. The disc sits at the end of the vaginal canal, directly below the cervix. It is therefore very important to push it in as far as possible and fix it in place. To keep it in place, position the outer edge of the disc, the one with the ring, behind the bone at the front of your pelvis, called the 'pubic bone'.

Once it's in place, it sits so comfortably in your body that you won't feel it at all. When positioned correctly, the disc sits at a diagonal angle inwards so that it sits in the right place for optimum function and comfort.

Unlike menstrual cups, it does not rely on the creation of a vacuum to stay in place, making it much more user-friendly and easier to use.

HOW TO USE IT?

Scan the QR code to watch the video guide on how to use, insert, and clean your menstrual disc.


The Femi.Eko® menstrual disc should be sterilised at the time of purchase before you use it for the first time. You can sterilise the disc by boiling it in water for 3 to 5 minutes (be careful

not to let the disc touch the bottom of the boiling pot) or by using a microwave, UV or other steriliser from www.femieko.com.


Hygiene is always important when we are talking about the genital area, and not just there! So wash your hands thoroughly before using the disc. Also wash the menstrual disc with pH-neutral, oil-free soap and water. Rinse thoroughly!


Fold the Femi.Eko® menstrual disc in the shape of an 8 and position the extraction ring towards the outside of the vagina. See? It's no bigger than a tampon, so relax, it'll be a :) If you wish, you can use water or a water-based lubricant (not silicone) to make insertion smoother.

Hold the disc firmly between your thumb and forefinger and start pushing from close to close through the opening of the vulva. As soon as you have passed the entrance to the vagina, push the disc further (still holding it firmly) and tilt it towards the coccyx. Push as deeply as possible until you feel that the disc can no longer penetrate anywhere.

Once you have pushed it in very deeply, it is important to position it over the pubic bone. Use your index finger to lift the edge of the ring disc up towards your navel. And you can forget about your menstruation for 12 hours (many women even report fewer cramps after using the disc).

To remove the disc, you must be relaxed. So don't forget to breathe :) The disc has nowhere to go, so you can get it out! Insert your finger into your vagina and grasp the special ring of the disc. Grasp the ring and gently pull the disc out. If you can't reach the disc, read the tips on how to use it.

Empty the disc (we recommend doing this directly after removing it into the toilet by simply grabbing the ring without the risk of dropping it into the toilet), rinse it thoroughly with cold water and wash it with neutral soap,

rinse again. Then you're ready to use it again. longer penetrate anywhere. When the cycle is finished, wash the disc, sterilize it, and store it in the special Femi.Eko® bag you received in the package.

TIPS FOR USE:

- Empty it at least twice a day (every 12 hours); we recommend before going to bed and after waking up.
- Always wash your hands and the menstrual disc thoroughly before and after handling.
- Do not share the menstrual disc with other people.
- Replace the menstrual disc if it is torn or damaged.
- The menstrual disc is not recommended for use during bleeding after childbirth.
- You can use the menstrual disc with an intrauterine device or a contraceptive ring, but ask your doctor first.

- If you drop the menstrual disc in the toilet, you do not need to replace it. Before using it again, disinfect it by boiling it in water for 4-5 minutes (but no longer than 7 minutes) or wipe it thoroughly with a soft cloth and 70% isopropyl alcohol and allow the alcohol to evaporate.
- The menstrual disc does not protect against pregnancy or sexually transmitted infections (STIs). It is safe to use with most contraceptives, but ask your doctor if you have any questions.
- Wash the menstrual disc after sex.
- If the disc moves upwards during arousal or intercourse, allow it to return to a lower position before removing it.
- If you feel uncomfortable during intercourse, we recommend that you change the position of the menstrual disc or try a different position that is more comfortable for you and your partner.
- If you are unable to remove the disc, sit in a seated, squatting position and push with your pelvic muscles, as if you were having a bowel movement. This movement shifts its position and you can grasp the disc by the ring. If you cannot grip the ring, rotate the disc with your finger on the inside of your vagina until it reaches a comfortable position. If you cannot grip the ring, grasp the edge of the disc and pull gently.
- Do not injure yourself, be gentle with the movements you make and do not insert any objects that could injure your mucous membrane.

If you notice any leaks when using the disc, check whether you have inserted the disc deep enough or whether you have an anatomical feature, such as a retroverted uterus, that requires special adaptation. Make sure that the edge of the disc is positioned correctly over the pubic bone with the ring facing upwards to prevent slippage (the disc may not sit properly and slip). If you experience this repeatedly, you may want to change it to a bigger size.

- If you leak into the toilet when urinating or defecating, do not worry. This phenomenon is called -emptying and is a feature appreciated by many women. Adjust the disc behind the pubic bone afterwards, but remember to replace it at the recommended 12-hour intervals. This phenomenon does not occur when walking or during everyday activities, but is a natural process caused by the specific contractions that occur during urination.
- We recommend removing the disc over the toilet or in the shower to avoid splashing blood.
- If you show signs of menstrual anxiety- menophobia, we recommend that you use the Femi.Eko universal applicator to insert and remove the menstrual disc to make the experience easier and reduce contact with menstrual fluid.

CARE:

- Wash thoroughly after each use with warm water and a mild, pH-balanced intimate soap that does not contain silicone, oils or chemical ingredients. Avoid strong detergents and anything that could irritate the skin. Rinse thoroughly with warm water!
- For a vegan and natural alternative, you can mix a small amount of baking soda or apple cider vinegar with water, and use it to clean the disc, and then rinse it thoroughly with water before reusing it.
- Use a soft brush for dried blood deposits.
- Wash and iron the storage bag occasionally to ensure complete hygiene.
- Dry the menstrual disc thoroughly and store it in the cotton bag (or other permeable container). Never use a plastic bag or airtight container for storage.
- Some discolouration may occur over time. This is normal and may vary from person to person.

To avoid this, always rinse the cup with cold water first.

- If the menstrual disc cracks or shows signs of damage to the silicone (which can occur after using non-recommended cleaning products), it is time to replace the disc.
- Do not wash the menstrual disc with substances that can irritate the vaginal area. Avoid using cleansers, soaps with strong fragrances or colourings that can disrupt the natural pH or damage the silicone, such as silicone-based products, oils, etc.
- Do not soak the menstrual disc in liquids other than water. Silicone is a porous material and can absorb oils or detergents if used over a long period of time. These substances can irritate the vaginal area or damage the silicone.

SAFETY AND PRECAUTIONS:

Keep the menstrual disc away from children and pets.

Never use the menstrual disc for more than 12 hours at a time.

If you experience symptoms such as general burning, irritation or swelling in the genital area or discomfort when urinating, you should stop using the disc and consult a doctor.

The menstrual disc does not protect against pregnancy or sexually transmitted diseases.

Internal menstrual products have been associated with Toxic Shock Syndrome (TSS). TSS is a rare but serious condition that can progress rapidly and lead to death. If you experience fever, severe abdominal pain, rash or vomiting while using the disc, stop use immediately and seek medical advice. By reading this information, you are ensuring that you are aware of the signs of TSS and if anything unusual happens, contact your doctor as soon as possible.

QUESTIONS?

If you have any further questions or concerns, please do not hesitate to contact us. We are here to support you and provide you with the information you need to have the most enjoyable and healthy experience possible. Contact us at: hello@femieko.com

We will be happy to answer any questions you may have about the menstrual disc.

Thank you for choosing Femi.Eko® - for a healthier and more environmentally friendly menstrual experience!

Manufactured in P.R.C., for FEMIEKO APS, Ryttervangen 11, 7323 Give, Denmark

Distributed with love by FEMIEKO GROUP S.R.L., DV24, C1-24, Bucharest, Romania

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DACH Distribution: Abraxa Business Network GmbH, Am Schafberg 2a. 96489 Niederfüllbach, Deutschland.