pepita

P-CMC

COMFORT Massage Chair –Instructions for use

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








Step 2

Step 4

Thank you for choosing our product! Please read this manual carefully before installing and using the armchair to ensure a safe and maximally relaxing experience.

Health and Usage Restrictions

To protect your health and safety, please consult your physician before using the massage and heating functions if you have any of the following conditions:

- **Cardiovascular diseases:** Especially an implanted pacemaker, severe arteriosclerosis, or deep vein thrombosis.
- **Pregnancy:** The use of the vibration function (especially in the lumbar area) is not recommended during pregnancy.
- **Spinal disorders:** In case of acute herniated disc, osteoporosis, or recent spinal injury.
- **Malignant tumors:** Vibration that stimulates microcirculation should be avoided in case of cancerous diseases.
- **Sensory loss:** If the user has reduced heat sensation (e.g., neuropathy due to advanced diabetes), the use of the heating pad is prohibited or requires extreme caution to avoid burns.
- **Acute inflammation:** Massage is not allowed in the presence of fever, acute inflammatory processes, or fresh surgical wounds.
- **Age restriction:** The device must not be used by children. Do not allow them to play with the remote control or the reclining mechanism.
- **Persons with reduced capabilities:** Persons with limited physical, sensory, or mental capabilities may only use the chair under supervision.
- **Weight limit:** The manual mechanism and frame structure of the chair can support a maximum body weight of up to 150 kg.
- **Use while sleeping:** Do not use the massage function while sleeping. Always set the TIMER function if there is a risk of falling asleep.

Technical Data and Specifications

- **Power supply (Input):** AC 100-240V, 50/60Hz
- **Power supply (Output):** DC 12V, 1A
- **Compliance standards:**
EN IEC 55014-1:2021,
EN IEC 61000-3-2:2019+A1:2021,
EN 61000-3-3:2013+A2:2021+AC:2022-01,
EN IEC 55014-2:2021
- **Product dimensions:** 88x91x102 cm / 34.1 kg
- **Max. load capacity:** 150 kg
- **Massage system:** 4-zone vibration motor (Back, Lumbar, Thigh, Leg/Calf)
- **Programs:** 5 different massage modes
- **Intensity:** 2 levels (Low / High)
- **Timer:** 3-level automatic shut-off (15 / 30 / 60 minutes)
- **Heating:** 1 separate, independently switchable heating pad

Assemble the massage chair based on the illustrated assembly instructions

Installation and Connection

- **Choosing the location:** Ensure sufficient space behind and in front of the chair for the reclining backrest and the footrest.
- **Power supply:** Connect the power adapter to the cable coming out of the chair, then plug it into a standard wall outlet.
- **Remote control:** Locate the wired remote control placed in the pocket.

Functions of the Remote Control Buttons

Using the wired remote control, you can precisely adjust the massage and heating functions of the chair:

1. **POWER:** Start and stop the massage function.
2. **INTENSITY:** Adjust the vibration strength in two levels (HIGH / LOW).
3. **HEAT:** Turn the built-in heating pad on and off separately.
4. **MODE:** Switch between the 5 different massage programs:
 - **Pulse:** Pulsing massage
 - **Press:** Pressing-type massage
 - **Wave:** Waving vibration
 - **Auto:** Automatically changing program
 - **Normal:** Traditional, continuous vibration
5. **BACK:** Switch for massaging the back zone.
6. **LUMBAR:** Switch for massaging the lumbar zone.
7. **THIGH:** Switch for massaging the thigh zone.
8. **LEG:** Switch for massaging the leg/calf zone.
9. **TIMER:** Set the automatic shut-off time (15, 30, or 60 minutes).

Operating the Manual Reclining Mechanism

Since the chair has a mechanical structure, you can adjust the recline angle using your own body weight and the side lever:

- **Opening the footrest:** Pull the lever on the side of the chair to release the footrest.
- **Reclining the backrest:** After the footrest pops up, shift your center of gravity backward and gently push the backrest back to the desired position.
- **Resetting to default position:** Lean forward so that the backrest straightens up, then firmly push the footrest back under the chair with your foot until it clicks into place.

Safety Rules and Operation

- **Time limit:** To prevent the motors from overheating, use the massage for a maximum of 30 minutes at a time. Allow the system to rest for at least 15 minutes before restarting.
- **Getting out:** Never try to stand up from the chair with the footrest open or the backrest reclined. Always return the chair to its default position before getting out.
- **Environment:** Ensure at least 50-60 cm of free space behind the chair to the wall so that the backrest can recline without obstacles. Use it on flat, even ground.
- **Children and pets:** When closing and opening the mechanical scissor structure, make sure that small children or pets are not in the immediate vicinity of the chair.
- **Power supply:** Use only the original factory power adapter included with the product.

EN

Maintenance and Cleaning

- **Disconnecting power:** Always unplug the power cord from the wall outlet before cleaning, moving, or maintenance.
- **Moisture:** Wipe the upholstery only with a slightly damp, thoroughly wrung-out soft cloth. Make sure no water gets into the electrical parts.
- **Chemicals:** Do not use harsh abrasives, bleach, or solvents.
- **Vacuuming:** Remove dust accumulated along the seams using a vacuum cleaner set to low suction power with a soft brush head.
- **Cable protection:** Regularly check that the remote control cord and the power cable run freely and are not pinched under the metal scissor structure.

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