
pepita

P-M67PRO-FBMC

M67 PRO Massage Chair –Instructions for use

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DE

Introducing the Pepita M67 PRO Massage Chair, delivering a revolutionary new massage experience right in the comfort of your home.

Crafted from high-quality materials, this chair provides professional-level full-body massage, while its customizable programs ensure that everyone finds their perfect massage mode. The premium PU leather upholstery and thick, plush padding guarantee ultimate comfort.

With the Zero Gravity function, the back and spine are completely relieved of pressure, allowing the entire body to relax deeply. By reducing the gravitational load on the torso and vertebrae, it alleviates back pain, enhances the effectiveness of the massage, and improves blood circulation. Ultimate relaxation is ensured by the 8D fixed-point massage system, soothing airbags, and a built-in heating function.

Enjoy peaceful pampering and the beneficial, healing effects of a professional massage in the embrace of the Pepita M67 PRO Massage Chair!

**Please read this manual carefully before operating the product.
Keep this manual for future reference!**

Important Safety Warnings

Before using this product, please read all instructions first.

To reduce the risk of burns, fire, electric shock, or personal injury, follow these instructions:

1. Do not use a power supply exceeding 220V AC.
2. Always unplug the power cord from the socket when the device is not in use daily, or before disassembling/repairing any parts.
3. To avoid fire, electric shock, or injury due to overheating, do not operate the product under pillows, blankets, or other coverings.
4. When the product is used by or near children, persons with physical or mental disabilities, or the disabled, supervision by a responsible adult is required.
5. Only use this product for the purposes described in the user manual. Do not use accessories not recommended by the manufacturer.
6. Do not use the product if any of the following occur: the power cord or plug is damaged, the unit malfunctions, it has been dropped or damaged, or it has been immersed in water. Return the product to an authorized service center for inspection and repair.
7. Do not pull the power cord to move the product, or use the power cord as a handle.
8. Keep the power cord at a safe distance from heated surfaces.
9. Do not operate the product with blocked air vents. Keep vents clear of lint, hair, and other debris.
10. Do not insert or drop any objects into any openings of the product.
11. Do not use the product outdoors.
12. Do not use the product in areas where aerosol products are being used or where oxygen is being administered.
13. To disconnect the power, turn off all controls first, then unplug the power cord from the socket.
14. Always supervise and ensure children do not treat the massage chair as a toy, to ensure their safety.

EN

Safety Instructions

1. Safety Precautions

- Do not use the massage chair in excessively humid or dusty environments, as this may cause electric shock or equipment malfunction.
- Do not use the massage chair where the room temperature exceeds 40°C.

- Do not place the massage chair near heaters, stoves, or in direct sunlight.
- This massage chair is designed for personal use only and is not suitable for any commercial purposes.

Caution

- The massage chair is for indoor use only.
- When using the massage chair in a low-temperature room, do not increase the temperature suddenly. It is recommended to gradually adjust the temperature to a normal range before use.
- If the massage chair is moved from a low-temperature storage area to a warm environment, moisture may accumulate inside the unit, which could affect its normal operation and even cause equipment failure.

Caution

- This product is not intended for self-diagnosis or treatment. If diagnosis or treatment is required, it should be handled by professional medical personnel.
- **Do not use this product if you:** have a cardiac pacemaker or other implanted medical electronic devices, are pregnant, have recently given birth, or are menstruating, suffer from high blood pressure, heart disease, cerebral arteriosclerosis, or cerebrovascular disease, have a malignant tumor.
- If you are currently taking medication or have a medical condition, consult your physician and follow their advice before using this product.

2. Safety Precautions

- To ensure safe and correct use of the massage chair, please read the operation manual carefully before use.
- To avoid reduced massage effectiveness or injury, do not use the massage chair simultaneously with electric blankets or other therapeutic devices.
- Do not use the massage chair within 1 hour after eating, after drinking alcohol, or when feeling unwell.
- Do not use the massage chair with wet hands or other wet parts of the body.
- Recommended usage: Maximum twice daily, 15-30 minutes each session. When not in use, turn off the power switch or unplug the chair. Before reuse, ensure the chair has been unplugged and inactive for at least 10 minutes. Prolonged continuous use may trigger the overheat protection device, causing an automatic power cut.
- If any liquid is spilled or poured onto the massage chair, immediately unplug the power cord from the socket.
- Do not operate the massage chair if the power cord, power switch, or power switch housing is damaged, as this may cause electric shock or equipment malfunction.
- Do not operate the massage chair if it has been damaged by impact during movement or use, or if water has entered the unit.


Situations Where the Massage Chair Should Not Be Used

If any of the following occurs, turn off power, unplug, and contact customer service immediately.

- If any liquid is spilled or poured onto the massage chair, as this may cause electric shock.
- If the upholstery of the massage chair is damaged, exposing internal components.
- If you experience any discomfort or severe pain during massage, stop the operation immediately and seek medical advice.
- If you notice any equipment malfunction or abnormal conditions during use, stop using the chair immediately.

- If there is a power interruption during use, take precautions to avoid injury when power is restored suddenly.
- During thunderstorms.

Part Diagram

Product Specifications

- Product Type: M67 PRO
- Rated Voltage: 220-240 V
- Rated Frequency: 50 Hz
- Rated Power: 90W
- Safety Class: Class I
- Product Dimensions: 142x75x110 cm
- Product Weight: 53 kg
- Material: ABS plastic, PU leather
- Maximum Weight Capacity: 200 kg

Compliance Standards:

EN 62479:2010

EN 50663:2017

EN 60335-2-2:2003+A1:2008+A2:2015

EN 60335-1:2012+A11:2014+A13:2017+A1:2019+A14:2019+A2:2019+A15:2021

+A16:2023

ETSI EN 301489-1 V2.2.3 (2019-11)

ETSI EN 301489-17 V3.3.1 (2024-09)

ETSI EN 300328 V2.2.2 (2019-07)

Features and Functions

- 8D fixed-point massage system (neck, back, waist)
- 5 automatic massage programs
- Automatic body scanning technology
- AI intelligent voice control
- Full-body airbag massage (arms, shoulders, calves, feet)
- Foot roller massage
- Heating function
- Bluetooth connectivity
- Touchscreen control panel
- Zero Gravity function
- Airbag intensity – 3 different levels
- Massage speed – 3 different levels
- Time setting: 15 – 30 minutes
- Decorative capsule design
- Easy and convenient operation
- Quiet operation
- Skin-friendly PU leather upholstery
- Convenient movement with casters


Pre-Use Instructions

1. Unpacking & Setup

This product is delivered as a single integrated unit. After unpacking, carefully move the massage chair into position. No assembly is required – simply connect the power cord to start using.


2. Pre-Use Preparation

After removing the packaging, locate the main power switch on the rear of the massage chair. Connect the power cord, then turn on the switch to complete the setup and prepare the chair for use.


3. Placement of the Massage Chair

Do not place the power cord under the massage chair's control box. (See Diagram 01/02)

Correct (Diagram 01)


Incorrect (Diagram 02)

1. Ensure adequate space around the massage chair before use (See Diagram 03). The installation area must provide sufficient space for the backrest to recline and the footrest to extend. Maintain a


minimum clearance of 35 cm behind the chair and 35 cm in front of it. Ensure no obstacles or people are in front of or behind the massage chair. To prevent floor damage, place a rug under the massage chair.

4. Moving the Massage Chair

Press both hands against the backrest to tilt the chair backward by approximately 45°, then carefully push or pull it forward or backward using the casters until it is in the desired position. (Diagram 04)


(Diagram 03)


(Diagram 04)

Control Panel and Operation

- **Power button:** Press to activate when powered; press again to disable all functions.
- **Shoulder & Hip airbag:** Press this key to turn on/off shoulder and hip airbag.
- **Hand airbag:** Press this key to turn on/off hand airbag.
- **Leg airbag:** Press this button to turn the leg airbag massage on/off.
- **Air pressure +:** Press "Air Pressure +" to raise airbag intensity (airbag must be on).
- **Air pressure -:** Press "Air Pressure -" to lower airbag intensity (airbag must be on).
- **Kneading mode:** Press this button to turn the neck kneading mode on/off.
- **Tuina mode:** Press this button to turn the back acupressure mode on/off.
- **Forward reverse:** Press this button to turn the forward and reverse massage function on/off.
- **Speed +:** In manual mode, press "Speed +" to increase kneading/acupressure intensity (requires mode to be on).
- **Speed -:** In manual mode, press "Speed -" to decrease kneading/acupressure intensity (requires mode to be on).
- **Fully automatic:** There are 5 automatic modes in total
 - F1. Relax Mode
 - F2. Strong Mode
 - F3. Shoulder & Neck Mode
 - F4. Waist & Back Mode
 - F5. Sleep Mode
- **Timing:** Press to set massage time (default 15 min, adjustable 5–30 min).
- **Pick up:** Hold "Pick Up" to raise backrest; release to lock angle.
- **Tilt down:** Hold "Tilt Down" to lower backrest; release to lock angle.
- **Hip massage:** Press to cycle hip roller intensity: Level 1 → Level 2 → Level 3 → Off.
- **Heating:** Press this button to turn the heating function on or off.
- **Foot roller:** Press to cycle foot roller intensity: Level 1 → Level 2 → Level 3 → Off.
- **Zero gravity:** Press this button to change the body angle and enter the zero-gravity posture. Press again to exit zero-gravity.
- **Airbag automatic:** Press this button to turn the full-body airbag massage function on/off.
- **AI voice button:** Press this button to turn the AI function on/off.

Shortcut Keys and Operation

Button	Description
	Power Key: Long press to turn the massage chair on or off.
	Knob: Rotate to adjust the massage intensity (Level 1, 2, or 3)
	Backrest Lift Key: Long press to raise the backrest.
	Backrest Recline Key: Long press to lower the backrest.
	Zero Gravity Key: Press to turn the Zero Gravity mode on or off.
	Auto Mode: Press to switch between different automatic modes.
	Airbag Key: Press to turn the full-body airbag massage on or off.
	Pause Key: Press to pause the massage, and press again to resume.
	Manual Mode: Press to switch between different manual modes.
	Heating Key: Press to turn the heating function on or off.
	Foot Roller: Press repeatedly to turn on the foot rollers, adjust the intensity (Level 1, 2, 3), or turn them off.
	USB charging port

AI Voice Control Function

No.	Command word
1	HI-MASSAGE-CHAIR
2	START-MASSAGE
3	END-MASSAGE
4	INCREASE-MASSAGE
5	REDUCE-MASSAGE
6	START-MODE-ONE
7	START-MODE-TWO
8	START-MODE-THREE
9	START-MODE-FOUR
10	START-MODE-FIVE
11	START-MODE-SIX
12	START-HEAT

No.	Command word
13	END-HEAT
14	INCREASE-AIR
15	REDUCE-AIR
16	START-FOOT
17	END-FOOT
18	BACK-RISE
19	BACK-DOWN
20	BACK-STOP
21	START-AIR
22	END-AIR
23	START-ZERO
24	END-ZERO

Bluetooth Music Playback

- Turn on the Bluetooth function on your mobile phone.
- Find the Bluetooth name of the armchair: M67PRO
- Enter the password: 8888, then use the music player to play music.

Daily Maintenance

1. For daily cleaning, use a neutral detergent and a soft cloth. Do not use corrosive detergents.
2. If the power cord is damaged, it must be replaced by a professional from the manufacturer's authorized service center to avoid hazards.
3. Do not use sharp objects to scratch or cut the surface of the product.
4. Unplug the power cord when the product is not in use or when performing cleaning to avoid hazards.
5. Do not place the product in a humid environment. Keep it in a dry and well-ventilated place.
6. Do not apply excessive force or heavy pressure to the product.

Troubleshooting and Maintenance

<i>Failure phenomenon</i>	<i>Possible cause</i>	<i>Processing method</i>
Noise <ul style="list-style-type: none">• Operating noise• Airbag inflation sound• Airbag deflation sound• Solenoid valve switching sound• Air pump operating sound	Occurs when the product emits mechanical operating sounds.	This sound is normal and requires no action.
Does not work at all	<ul style="list-style-type: none">• The power plug is not inserted into the outlet.• The switch on the socket is not turned on.• The Massage button is not pressed.	Check the power plug and outlet, ensure proper contact, then press the power switch and the massage button.
The ongoing massage suddenly stops	The product has just finished its 15-minute rated working time and automatically switched to standby mode.	If you want another massage, press the ON button to activate the massage function again.
	The rated operating time has been exceeded, and the motor overheating protection has been activated.	Turn off the power, disconnect it from the power source, and wait for the motor to cool down before further use..

Environmental Protection Marking (10 years)

Under the normal operating conditions described in this manual, the environmental protection service life of the product is 10 years. It can be used safely within this period. After the expiration of the environmental protection service life, the product should enter the recycling system.

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