SmileHOME

by Pepita

SH-0027

PRO 300B FOOT AND CALF MASSAGER

EN

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HU

PRO 300B MASÁŽNY PRÍSTROJ NA NOHY A LÝTKA

SK

PRO 300B MASER PICIOARE SI GAMBE

RO

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DE

Thank you for choosing our product!

SmileHOME products are made with the highest quality and state-of-the-art technology.

Please read the following information carefully before use and keep the guide for future reference. The safe operation of the appliance can only be guaranteed if you follow the warnings and safety instructions concerning general safety and accident prevention.

No liability is accepted for damage resulting from improper use or operation!

IMPORTANT SAFETY REGULATIONS

WARNING

- Unplug the appliance after use or before cleaning the product.
- Do not cover it during use to avoid malfunction or damage due to overheating of the motor.
- Children, physically and mentally handicapped persons should use the appliance only under supervision.
- If you are under medical treatment, please consult your doctor before use.
- Please stop using the product immediately if you feel uncomfortable.
- If the power cords or plugs are damaged, the product is not working properly, or the product has been dropped, damaged or dropped into water, do not use it. Please contact your service centre for inspection and repair!
- Do not use spare parts not recommended by the manufacturer.
- Do not handle the button with wet hands.
- Keep the buttons away from water or other liquids.

OPERATING ENVIRONMENT

- Do not use the product in a wet or dusty environment to avoid product failure.
- Do not use the product at room temperature above 40°C.
- Do not use the product outdoors.
- Keep away from heat.

WARRANTY

The warranty does not cover the following situations:

- Damage caused by human factor, physical damage.
- Failure to follow instructions.
- Improper storage.
- Damage resulting from improper use.
- Dismantling, repair or modification of the product without the manufacturer's permission.
- Damage caused by users after purchase due to improper shipping.
- Damage resulting from any other force majeure (such as flood, fire, lighting and earthquake).

RELIEVING PAIN IN THE CALF

- Set the massager in the position shown in the picture.
- Place your calf on the massager.
- In this position the massager stimulates blood circulation and relaxes the calf muscle.

MASSAGING THE FEET

- Do not wear socks for this function.
- Move the foot slowly and rhythmically towards the arch.
- The increased blood flow will relax your legs.

You can use the massager in several optional positions:

HOW THE PRODUCT WORKS, REMOTE CONTROL OPERATION

Automatic program: P0-P9, Manual mode HR.


In manual mode you can adjust the intensity, in automatic mode you cannot.

Speed intensity levels: 1-30

Timing: 1-15 minutes

Power button

Change the speed
and time of the program.

time setting

manual mode

power button

speed setting

remote control sensor

TECHNICAL DATA

- Power: 48 W
- Voltage: DC 24 V
- Rated input: AC 100-240 V ~50/60 Hz
- Rated output: 24V 2A
- Material: ABS
- Display: LCD screen
- Timer: 1-15 minutes
- Speed intensity levels: 1-30
- Massage programs: P0-P9
- Product size: 37.5*31*25.5 cm
- It has a stress-reducing and pain-relieving effect.
- It relaxes muscles, soothes tired joints and improves circulation.
- With AC adapter and remote control (AAA battery. Not included.)

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