SmileHOME

by Pepita

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SMART PERSONAL SCALE

User's Guide

Please read the user manual carefully before the first use and keep it for future use.



IMPORTANT! The mobile phone application is not yet available in Hungarian, but it can be changed to the most common foreign languages (English, German, French, Spanish, Italian, etc.).

Since the products (scale and application) are constantly being updated and improved, the images in the user guide are for informational purposes only.


Specification

	Parameters
Name of Product	Smile Home Smart personal scale
Size	26 x 26 x 2,3 cm
Connection method	Wireless connection with mobile phone application
Power supply	USB charging
Display	LCD / LED
Unit of measure	Kg / LB
Measuring range	0.2 -180 kg
Divisor value	10 g
Low voltage sign	LO
An overload signal	OL

Product structure diagrams

Note: USB charging requires a data cable. Charging time is less than 2 hours.

Using the scale

Precautions before use:



1. Please place the scale on a flat and firm ground first. Avoid on soft ground use (for example: wooden floor, grass, carpet, etc.)!
2. Check that the power supply of the scale is correct! To do this with one leg put a little weight on the scale!
3. When the display shows “0.00”, the weight is deleted. Then when the screen goes black, will stop automatically, thus indicating the end of the operation.

Download and install the application on your Android / iOS phone!

Method 1: Scan the QR code attached to the product for installation!

Method 2: Use your mobile web browser or Google Play / App Store download the OKOK International application from the online store!

You can register and log in to the application with your e-mail address and phone number.


If you cannot log in to the application:

1. Make sure that your phone's internet connection is unobstructed!
2. Please contact the manufacturer to check if the app handles it the updates!
3. It is possible that the username has already been registered by someone else. (If you have excluded the above two points, we recommend that you enter a number after the registered name!)


Click on the Add Device icon!


Connect your phone to the scale via Bluetooth („bind/pair”)!


Turn on the „Bluetooth” and „Location Services” services on your mobile phone!

Select the device and connect it!

Look for it on the display, then press the scale icon!


Weight measurement

1. When measuring, turn on the „Bluetooth“ connection and the „Location“ service on your phone!
2. Remove the protective film from the surface of the 4 electrode sheets.
3. Then press the scale icon (weigh) again!
4. Please stand on the scale!

- last time
- best in 30 days
- target weight
- history
- weight
- BMI
- fat
- muscle
- water
- visceral fat
- bone mass
- metabolism
- obesity degree
- weight without fat

The appearance of this interface indicates that the test was successful. Click on it and you can see the measurement data („bind/pair“)! If the test is unsuccessful, please check whether you have followed the prescribed usage methods. (Dry skin can also affect the measurement. You can try the measurement again by washing and then semi-drying your soles!)

Maintenance

1. The surface of the scale can be cleaned with a slightly damp cloth. Do not use corrosive cleaning products!
2. Please do not put the product in water!
3. To avoid scratches, please do not place the scale with the display surface facing down!
4. Please keep the device away from babies and small children.
5. The product is made of glass, so be extra careful! Protects it from impacts and possible-falling!
6. Store in a dry, well-ventilated, cool place!
7. The product is intended for home use and cannot be used for commercial measurement.

Frequently Asked Questions

1. Why is the later measurement result different from the initial one?

After the first use or after moving, the first measurement data also contain the weight of the product. Please wait until the device switches off automatically before taking another measurement!

2. Why is there such a big difference in the measurements on different soils between?

Due to the sensitivity of the electronic scale, be sure to place it on solid and flat ground when using it! If the ground is uneven or soft objects are placed on it (such as a carpet), the measurement data will be inaccurate.

3. Why are the measurement data different at different times on the same day?

It is normal if the body weight is different in the morning and in the evening. Scientific observations show that due to metabolism and other factors, people's body weight can change by about 3.5 kg in a day. Changing clothes, changes in the body before and after meals, and before and after training are associated with changes in weight.

In what cases can the measurement results change?

Before and
after training

Morning and
evening error

Posture

Soil factor

Before and after
going to the
bathroom

Increasing and
decreasing
garments

Before and after meal

Loading

Safety regulations

ATTENTION! The scale is made of tempered glass. For the safety of you and your family's health and property, please take the following into account!

1. Do not heat used batteries or throw them into fire to avoid explosion!
2. We ask customers with reduced mobility to take the measurements with the help of others or use auxiliary devices (e.g. handrails)!
3. Do not measure in the dark or on wet ground!
4. Do not place objects under the scale!
5. Please perform the measurement on flat and hard ground!
6. To avoid accidents, avoid jumping when measuring!

1. Pregnant women and the elderly over 70 years old, as well as patients treated with fluoroscopy and patients with symptoms of edema may have different body fat ratios.
2. In the case of minors, due to their physical characteristics, the measurement will probably be inaccurate, because there is not an adequate amount of data available.
3. The use of the product is not recommended for people with pacemakers or other implanted medical devices.

Guarantee

Complaint/complaint

Complaints and complaints can be submitted on the www.pepita.hu website, cancellation and guarantee, or by filling out the appropriate form under the complaint menu items.

Guarantee

The product is covered by a 1-year manufacturer's warranty if the product was used as intended and in accordance with the regulations.

Place of origin: China

Manufacturer: Guangdong ENNO Medical Technology Co., Ltd.

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Výrobca: Guangdong ENNO Medical Technology Co., Ltd.