SmileHOME
by Pepita

SH-MPLY-2-6SL

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EN

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HU

MINI MASÁŽNA PIŠTOĽ

SK

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RO

MINI-MASSAGEPISTOLE

DE

Thank you for choosing our product!

SmileHOME products are made with the highest quality and state-of-the-art technology.

- With 4 different head designs and 6 speed settings, you can massage your muscles to perfection.
- Its ergonomic design ensures a perfect massage experience on any part of the body.
- It reduces muscle pain in a short time, relaxes and soothes.
- Effectively relieve pain caused by lactic acid build-up after intense exercise, it has a positive effect on the circulatory system, reducing the pain of muscle injuries and strains.
- Its practical battery design makes it easy to transport and use.

Please read the following information carefully before use and keep the instructions for future reference.

Safe operation of the device can only be guaranteed if the general safety and accident prevention warnings and safety instructions are observed.

No liability is accepted for damage resulting from improper use or operation!

TECHNICAL DATA

- | | |
|---------------------|-----------------------|
| • Display: | LCD display |
| • Speed: | 6 levels |
| • RPM: | 1800-3200 rpm/min |
| • Product size: | 14x4.5x13.5 cm |
| • Product weight: | 500 g |
| • Voltage: | 5 V |
| • Battery capacity: | 1200 mAh |
| • Charging way: | USB-Type-C, 8.4 V/1 A |
| • Box size: | 19.5x19.5x6 cm |

BOX CONTENTS

- 1x mini massage gun
- 4x replaceable heads
- 1x USB-Type-C charging cable
- 1x carrying bag
- 1x user manual

EN

CHARGING

- For the first time, please charge the battery for 3 hours.
- Please connect the data cable to the 5V adapter for charging. When charging, the charge indicator displays the electric quantity, while working, it displays the speed gear.
- The battery can be recharged at any time and at any charge level. Do not use the device while charging!

OPERATION

- Power ON/OFF button on the display - long press to start the massager / long press to turn it off.
- Press the + or - key to change the speed level from level 1 to level 6.
Each press increases or decreases the level by one.
- Detachable and replaceable massage head. Switch off the machine before connecting /changing the massage head. When disassembling the massager, just pull out the head with force, and to connect, push it in until it stops.
- Please select 5V charging adapter for charging (not included).

4x replaceable heads

USING MASSAGE HEADS

SPHERE MESSAGE HEAD

Due to its large surface area, it enables targeted massaging of the entire body and larger muscle groups, such as arms, back, glutes, hamstrings, legs, etc. We also use this in case of surface sensitivity due to injury.

FORK MESSAGE HEAD

Designed for massaging long and narrow muscles, such as the muscles around the spine (be careful not to touch the edges of the vertebrae), the neck and the Achilles tendon.

CONE-SHAPED MESSAGE HEAD

It is suitable for deep tissue massage, for example, around the shoulder blades, joints, palms, legs, and soles. Avoid contact with bones!

SMALL, MUSHROOM-SHAPED MASSAGE HEAD

It is designed for smaller surfaces, so it is also perfect for massaging denser muscle groups.

MAINTENANCE, CLEANING, STORAGE

- Clean the appliance with a slightly damp cloth, then wipe dry with a soft cloth.
- Store in a dry, well-ventilated place free from corrosive gases.
- If the appliance is not used for a long period, recharge every six months.

WARNINGS AND PRECAUTIONS

- Do not place in water. Keep away from liquids or heat sources. Keep ventilation openings clear.
- Do not remove screws or attempt to disassemble.
- Do not use continuously for more than one hour. Allow the machine to rest for about 30 minutes before use again.
- Disconnect from the charger before use.
- For adult use only. Do not use if the device is damaged. In case of illness, consult your doctor before using the product.
- Use only on a dry, clean body surface, applying light pressure and for approximately 60 seconds per region.
- Apply only to soft tissue, according to complaints of pain or discomfort. Do not use on hard bony areas such as the skull.
- Keep hair and fingers away to avoid risk of pinching.
- Check the device carefully before each use.
- Never leave the device unattended during operation or charging.
- If a fault is detected, stop using the product and contact your service agent.

CONTRAINDICATIONS

- PLEASE DO NOT USE THE DEVICE OR OTHER VIBRATION DEVICES WITHOUT PRIOR APPROVAL FROM YOUR DOCTOR IF YOU HAVE ANY OF THE FOLLOWING: Pregnancy, diabetes with complications such as neuropathy (numbness) or retinal damage, while wearing a pacemaker, after recent surgery, epilepsy or migraine, herniated disc, spondylolisthesis, spondylolysis, spondylosis, endoprosthesis, intrauterine coil, metal screw or plate or any metal implant.
- Children and mentally challenged individuals can use this vibration device in the presence of an adult.
- The contraindications listed above do not exclude the possibility of using the device, but a consultation with a doctor is required before doing so. Research is currently ongoing in order to shorten the list of contraindications above. Practical experience shows that the integration of vibration treatment into the therapeutic plan is recommended in many cases.
- Before performing treatments, consult a doctor, specialist or physiotherapist.

ADVICE

The warranty does not cover the following situations:

- Damage caused by human factor, physical injury.
- Failure to follow instructions.
- Improper storage.
- Damage resulting from improper use.
- Dismantling, repair or modification of the product without the manufacturer's permission.
- Damage resulting from improper delivery by users after purchase.
- Damage resulting from any other force majeure (such as flood, fire, lighting and earthquake).

DISPOSAL INSTRUCTIONS

In order to protect the environment, the device must not be thrown away with household waste at the end of its useful life. Disposal must be carried out at the appropriate collection points in your country. Dispose of the device in accordance with the European Union Directive on Waste Electrical and Electronic Equipment (WEEE). If you have any additional questions, please contact the local organizations responsible for waste management.

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